

# Ultra-small wireless EEG registration modules, EEG-headsets and Bio-wristband

- Neurofeedback
- Biofeedback
- Neuromarketing

Ultra-small devices uses advance digital technologies for data acquisition and transmission (Bluetooth Low Energy) and expands perspectives for wide applications of distributed network registration of neurophysiologic and psycho physiological information about state of single person or group.

## EEG-registrators



## Bio-wristband module




for registration of photoplethysmogram, skin conductivity, electrocardiogram and movement activity



Elastic MCScap textile-caps with Ag/AgCl electrodes for contact gel (24 EEG channels)



# Wireless EEG recording using an EEG-headset and recording physiological signals using a Bio-wristband and PPG sensor



## Contact information

Frunze Str., 68, Taganrog, Russia, 347900  
**Phones:** +7 (8634) 62-62-42, 62-62-43, 62-62-44  
**e-mail:** office@medicom-mtd.com  
**www.medicom-mtd.com**